

# ONCOLOGY



ONCOLOGY

## HEPASPHERE™ MICROSPHERES

Doxorubicin-loaded HepaSpheres Microspheres have demonstrated success in the treatment of Hepatocellular Carcinoma (HCC) with high-objective tumour response while minimising the impact to underlying healthy liver function.<sup>1</sup> HepaSpheres Microspheres are reconstituted with saline or non-ionic contrast media and swell to approximately 4X their dry state diameter. These highly compressible microspheres facilitate outstanding contact with the vessel intima.




Catalogue Number	Size Range Dry Particles (µm)	Size Range When Expanded (µm)	Colour Code	Quantity (mg)
V225HS	30 - 60	120 - 240	Orange	25 mg
V325HS	50 - 100	200 - 400	Yellow	25 mg
V525HS	100 - 150	400 - 600	Blue	25 mg
V725HS	150 - 200	600 - 800	Red	25 mg

Packaged 1 per box.

HepaSpheres Microspheres are not approved or available for sale in the United States.

<sup>1</sup>Data on file.



*At the level of occlusion, the conforming nature of HepaSpheres Microspheres allows them to adapt to the tortuous anatomy of the vessel lumen providing permanent occlusion.*